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The screenshot shows the 'Functions and VIs' window with the 'GPIB Write' VI selected. The window title is 'Functions and VIs'. Below the title bar are three tabs: '索引(I)', '上一步(B)', and '列印(P)'. The main content area has a yellow background and contains the following text:

GPIB Write

Writes **data** to the GPIB device identified by **address string**.

Click the **parameters** for more information.

timeout ms (488.2 global)

The diagram shows the GPIB Write VI block with the following inputs and outputs:

- Inputs: timeout ms (488.2 global), address string, data, mode (0), error in.
- Outputs: status, error out.

The screenshot shows the 'Functions and VIs' window with the 'Concatenate Strings' VI selected. The window title is 'Functions and VIs'. Below the title bar are three tabs: '索引(I)', '上一步(B)', and '列印(P)'. The main content area has a yellow background and contains the following text:

Concatenate Strings

[String Function Descriptions](#)

Concatenates input strings and one-dimensional arrays of strings into a single, output string. For array inputs, this function concatenates each element of the array.

The diagram shows the Concatenate Strings VI block with the following inputs and output:

- Inputs: string 0, string 1, ..., string n-1.
- Output: concatenation of string0, string1, ..., string n-1.

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Functions and VIs

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
GPIB Read

Reads **byte count** number of bytes from the GPIB device at **address string**.

Click the parameters for more information.

timeout ms (400.2 global)

The diagram shows the GPIB Read VI block with the following connections:
Inputs: address string (blue), byte count (green), mode (U) (purple), error in (pink).
Outputs: data (blue), status (green), error out (pink).


Functions and VIs

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From Exponential/Fract/Eng

[String Function Descriptions](#)

Interprets the characters 0 through 9, plus, minus, e, E, and the decimal point (usually period) in **string** starting at **offset** as a floating-point number in engineering notation, or exponential or fractional format and returns it in **number**.

The diagram shows the From Exponential/Fract/Eng VI block with the following connections:
Inputs: string (blue), offset (green), default (0 dbl) (purple).
Outputs: offset past number (blue), number (green).

use system decimal point (T)
EFG

Note If you wire the characters Inf or NaN to string, this function returns the G values Inf and NaN, respectively.